Quality Specialty Products INC



Dormitory Chair HFC04 Specification









Features

- All materials be new and free of defects that may negatively affect appearance, durability, or function.
- 16 gauge x | .00 x 1.00" Square Tube Frame
- 16 gauge Back Mounting Brackets
- 16 gauge Seat Mounting Brackets

Frame

Each frame consists of a "H" shaped and "L" shaped piece of 16 gauge x 1.00" x 1.00" square tubing welded together to create the four legs. The seat and back is mechanically fastened with no less than eight (8) bolts for rigidity and strength. Plastic floor glides are attached to each leg to protect the chair and flooring. Plastic closure caps are inserted into the open ends of each tube.

Seat and Back Cushions

The seat consists of a 1/2" plywood seat support with 3" thick foam cushion. The back consists of a 1/2" plywood back support with 2.5" thick foam cushion. The seat and back are upholstered with fabric of sufficient quality to withstand high-use and hard-use environments

Quality Specialty Products INC













